


A.T.P. RING

CAPSULAR TENSION RING

Description

One piece PMMA ring
High flexibility due to special design
ETO sterilized



ATP CTRs are used to stabilize and recenter the capsular bag in the following cases:

- significant zonular weakness
- Lens subluxation
- Pseudoexfoliation
- Myopia
- Zonulolysis
- Soft IOL
- Marfan's Syndrome

Specification

Material: Blue PMMA/Clear PMMA

Model	Size	Thickness
CTR-12	10-12	0.20 mm
CTR-13	11-13	0.20 mm

